



Runs WatchBP Analyzer in Mac

WatchBP Analyzer are developed for Windows system.

There is an alternative solution For Mac user to run the WatchBP Analyzer by using Parallels Desktop for Mac.

<https://www.parallels.com/products/desktop/>

Fully licensed copy Parralels Desktop for Mac can be purchased online; alternatively you can download a free trial version.

www.parallels.com/en/buyonline

Microlife has tested the WatchBP Analyzer in Mac (Mac OS X El Capitan, version 10.11.2) using Parallels Desktop with Windows 10.

System Requirements

Hardware:

A Mac computer with an Intel Core 2 Duo, Core i3, Core i5, Core i7, or Xeon processor 4 GB of memory (8 GB recommended)

850 MB of disk space on the boot volume (Macintosh HD) for Parallels Desktop installation

Additional disk space for virtual machines (varies on operating system installed)

Internet connection

Software:

macOS Mojave 10.14 or later

OS X 10.11 El Capitan 10.11 or later

Disclaimer:

Microlife assumes no responsibility for errors or omissions in the use of Parallels software.